Published online 2019 September 28.

Brief Report

Colloidal Gold Assay, a New Screening Test for Detection of RSV Infection: A Controlled Trial

Lifang Liu¹, Wei Li¹, Ran Tao¹, Huamei Li¹ and Shiqiang Shang^{1,*}

¹Clinical Laboratory, Children's Hospital, Zhejiang University School of Medicine, Hangzhou, China

corresponding author: Clinical Laboratory, Children's Hospital, Zhejiang University School of Medicine, Hangzhou, China. Email: shangsq@zju.edu.cn

Received 2019 January 12; Revised 2019 May 26; Accepted 2019 May 31.

Abstract

Background: Respiratory syncytial virus (RVS) is the main cause of morbidity and mortality in the pediatric population. **Objectives:** The objective of the study was to compare the diagnostic accuracies of the colloidal-gold-Genesis and newly developed real-time RT-PCR-IngeniGen with an admitted real-time RT-PCR-HuaRuiAn as the gold standard for respiratory syncytial virus detection.

Methods: These three assays were performed on collected specimens.

Results: Altogether 258 specimens were studied. Among them, real-time RT-PCR-HuaRuiAn detected the respiratory syncytial virus in 166 samples, while colloidal-gold-Genesis detected 171 samples and real-time RT-PCR-IngeniGen detected 168 samples as positive. The overall sensitivity for the detection of the respiratory syncytial virus in comparison to real-time RT-PCR-HuaRuiAn was 95.2% for colloidal-gold-Genesis and 100% for real-time RT-PCR-IngeniGen. The specificity was 85.9% for colloidal-gold-Genesis and 97.8% for real-time RT-PCR-IngeniGen. The positive likelihood ratio was 6.75 for colloidal-gold-Genesis and 45.45 for real-time RT-PCR-IngeniGen. The negative likelihood ratio was 0.06 for colloidal-gold-Genesis and 0.00 for real-time RT-PCR-IngeniGen. The cycle threshold value of real-time RT-PCR-IngeniGen was higher than that of real-time RT-PCR-HuaRuiAn.

Conclusions: Colloidal-gold-Genesis is a good alternative for fast screening of respiratory syncytial virus detection.

Keywords: Respiratory Syncytial Virus (RSV), Real-Time RT-PCR, Colloidal Gold Assay

1. Background

Respiratory syncytial virus (RSV) is the most common cause of viral acute lower respiratory tract illness in pediatric population younger than 5 years (1-4).

Due to the identical clinical symptoms of RSV and other viral respiratory infections, an early, rapid and accurate diagnosis of RSV infection in susceptibile population is of great important (5, 6). On the other hand, the rapid diagnosis of RSV infection can exclude bacterial infection, minimizing inappropriate uses of antibiotics (7). There are a variety of approaches for detecting RSV infection, including antigen detection testing, cell culture, and PCR based methods worldwide (8-11). Cell culture and PCR are considered as the gold standard, but there are time-consuming, technically demanding, laborious, and expensive.

In Hangzhou, China, most hospitals are likely to screen RSV infection in outpatients by using colloidal gold assay which works by detecting the antigen of RSV and is characterized as quick detection, while the real-time RT-PCR assay usually functions as the confirming test combined with clinical symptoms.

2. Objectives

The objective of this study was to evaluate the diagnostic accuracies of the colloidal-gold-Genesis, (CFDA: 20143400267) and the newly developed real-time RT-PCR-IngeniGen in a head-to-head comparison with an admitted real-time RT-PCR-HuaRuiAn (CFDA: 20163400961) as the gold standard.

3. Methods

3.1. Study Population and Specimens

A prospective study was conducted from July 2017 to December 2017 in the Children's Hospital, Zhejiang Uni-

Copyright © 2019, Author(s). This is an open-access article distributed under the terms of the Creative Commons Attribution-NonCommercial 4.0 International License (http://creativecommons.org/licenses/by-nc/4.0/) which permits copy and redistribute the material just in noncommercial usages, provided the original work is properly cited.

versity School of Medicine. We collected the nasopharyngeal specimens regarding following conditions: (1) patients younger than 14 years old, who were suspected of having RSV infection with symptoms like cough, dyspnea, and wheezing (12), (2) the specimens were divided into RSV positive or negative group detected by colloidal gold assay, (3) then all three RSV detection methods were performed on the specimens by performers blinded to the former test results.

This study was approved by the medical ethics committee of the Children's Hospital of Zhejiang University School of Medicine (No. 2016-IEC-013, date: 2016.11.30) and written consent was obtained from parents before enrollment.

3.2. Colloidal-Gold- Genesis for RSV Detection

Colloidal gold plate (Genesis company, Hangzhou, China) contained three major parts: (1) sample well, (2) quality control part, (3) RSV test part. Specimens were tested according to the manufacturer's instruction. Briefly, 3 drops (about 100 μ L) of the sample were added into a sample well and incubated for 15 minutes. Samples with lines occurred in the quality control part and RSV test part were determined as RSV positive; samples with lines occurred in quality control part while the RSV test part showed nothing were determined as RSV negative. Samples showing nothing in quality control part indicated an invalid test, and a second test plate or other test assay was necessary till the result was either positive or negative.

3.3. Real-Time RT-PCR Testing

Total nucleic acid extraction kit by magnetic beads (IngeniGen, China) was utilized to extract the nucleic acid from all specimens. For the commercial RSV nucleic acid detection assay (HuaRuiAn Biology and IngeniGen, China), positive control, negative control, and quality control were included in each run. Real-time RT-PCR was performed in the ABI7500 system according to each brand manufacturer's instruction (Table 1) (HuaRuiAn Biology: http://www.huayinbio.net/product-4-92-158.html; Ingeni-Gen: http://www.ingenigenbiotech.com/col.jsp?id=106). Cycle threshold (CT) value of each sample was collected for future analysis.

3.4. Statistical Analysis

The results were analyzed by kappa test (SPSS software, version 24.0). Kappa value greater than 0.8 was considered as a very good agreement (13, 14). The average, median, minimum and maximum CT values were determined. The

comparison of CT value was analyzed by independent or paired *t*-test. A two-tailed P value of less than 0.05 was considered to be statistically significant. All P values were rounded up to 3 decimal places.

4. Results

During the study period, 258 children were enrolled. Among 258 specimens, RSV nucleic acid detection assay (real-time RT-PCR-HuaRuiAn) detected 166 positive samples and worked as the gold standard. 158 specimens were accurately identified as RSV positive by colloidalgold-Genesis, 13 specimens were incorrectly identified as RSV positive, and 87 specimens identified as negative by colloidal-gold-Genesis while 8 of them were identified as positive specimens by real-time RT-PCR-HuaRuiAn (Table 2). The overall sensitivity and specificity of colloidal-gold-Genesis for RSV detection in comparison to real-time RT-PCR-HuaRuiAn were found to be 95.2% and 85.9%, respectively. And a high statistical consistency was found between the two RSV detection assays (kappa value = 0.82). The positive and negative likelihood ratios were 6.75 and 0.06 respectively.

5. Discussion

BD VeritorTM System RSV and Quidel[®] Sofia[®] RSV FIA are the widely used lateral flow digital immunoassay (DIA) tests for RSV viral antigens throughout the world. Neena Kanwar et al. compared the diagnostic accuracy of the two DIA tests with the real-time RT-PCR assay results as a reference, showing a sensitivity of 71.15% for BD VeritorTM System and 80.77% for Quidel[®] Sofia[®] RSV FIA system, and the specificity of 100% for both systems (15). The sensitivity and specificity of the nCounter (NanoString Technologies), which is a multiplex digital method of RNA quantification, in comparison to real-time RT-PCR were found to be 74.3% and 98.4% for RSV-A, 77.6% and 97.8% for RSV-B, respectively (16). Compared with the results of those studies, the colloidal-gold-Genesis showed higher sensitivity and lower specificity with a higher Youden index. Thus, the colloidal-gold-Genesis is a good alternative for the clinical quick screening of RSV infection for its high specificity, sensitivity and time-saving. Furthermore, the median CT value for the specimens in the colloidal positive group was significantly lower than the negative group. The higher the CT value, the lower is the RSV titer, thus the result tested by colloidal-gold-Genesis is not credible enough, especially

Real Time RT-PCR-HuaRuiAn 1- 37°C for 25 min, 94°C for 2 min					Real Time RT-PCR- IngeniGen 45°C for 10 min, 94°C for 10 min				
3- 94°C for 5 sec,	55°C for 35 sec, 32 cy	cles							
ble 2. Compariso	n of Colloidal Gold A	ssay and Real-Time R	T-PCR-IngeniGen Resu	ults with Reference to	the Real-Time RT-P(CR-HuaRuiAn for R	SV Detectio	on	
	n of Colloidal Gold A True Positives	ssay and Real-Time R False Positives	T-PCR-IngeniGen Resu True Negatives	ults with Reference to False Negatives	the Real-Time RT-P(Sensitivity,%	CR-HuaRuiAn for R Specificity, %	SV Detectio	on LR-	Карр
ble 2. Compariso Test Assay Colloidal gold									Kapp 0.82

Abbreviations: LR+, positive likelihood ratio; LR-, negative likelihood ratio.

when the RSV titer is low, which needs another test for confirmation, that's why colloidal-gold-Genesis is a fast screening assay rather than a confirmation assay.

The RSV nucleic acid detection assay (real-time RT-PCR-IngeniGen) detected a total of 168 positive samples with 2 of them determined incorrectly, and 90 identified correctly negative samples. The sensitivity and specificity were found to be 100% and 97.8% respectively (Table 1). The positive and negative likelihood ratio was 45.45 and 0.00 respectively. A higher statistical consistency was found between the real-time RT-PCR-HuaRuiAn and real-time RT-PCR-IngeniGen (kappa value = 0.983).

As shown in Figure 1A, the CT value of colloidal gold assay positive group was significantly lower than that of the negative group (P \leq 0.000). The median PCR CT value for the positive group was 24.18 (range 17.83 - 28.8), and for the negative group 27.45 (22.32 - 30.14). To determine the difference between the two real-time RT-PCR assays, the CT value was detected and analyzed. The median CT value of the 166 RSV positive samples which were detected by real-time RT-PCR-IngeniGen was 30.76 (range 25.72 - 37.42), while median CT value of real-time RT-PCR-HuaRuiAn was significantly lower (24.21, range 17.83 - 30.14; $P \le 0.000$) (Figure 1B). Coming to the 2 incorrectly identified positive samples by real-time RT-PCR-IngeniGen, the CT value was 38.05 and 39.13 respectively. The 2 incorrectly identified positive samples by real-time RT-PCR-IngeniGen showed that the result of a positive sample may be inaccurate as the CT value

higher than 38, suggesting that the outcome criteria of CT value should be more precise.

In future, the detection kit for RSV infection will develop towards fast, automatic, high sensitivity and specificity, as well as point-of-care testing (POCT). The colloidalgold-Genesis is effective for early screening of RSV infection, while the newly developed RSV nucleic acid detection assay (real-time RT-PCR-IngeniGen) which is performed to confirm RSV infection should be more accurate in the definition of positive samples.

Footnotes

Authors' Contribution: Shiqiang Shang and Lifang Liu designed the study. Wei Li, Ran Tao and Huamei Li performed the experiments. Lifang Liu verified and analyzed the data, and wrote the manuscript. All the authors read and approved the final manuscript.

Conflict of Interests: All of the authors have no interests conflict.

Funding/Support: This study was funded by the National Nature Sciences Foundation of China (No. 81701535).



Figure 1. Boxplot illustrating CT value of specimens tested by the three assays: colloidal-gold-Genesis, real-time RT-PCR-IngeniGen, and real-time RT-PCR-HuaRuiAn. A, the comparison of real-time RT-PCR-HuaRuiAn CT value between the colloidal gold assay positive group and negative group. B, the comparison between the real-time RT-PCR-HuaRuiAn Biology CT value and the real-time RT-PCR-IngeniGen CT value for the 166 positive samples.

References

- Scheltema NM, Gentile A, Lucion F, Nokes DJ, Munywoki PK, Madhi SA, et al. Global respiratory syncytial virus-associated mortality in young children (RSV GOLD): A retrospective case series. *Lancet Glob Health*. 2017;5(10):e984–91. doi: 10.1016/S2214-109X(17)30344-3. [PubMed: 28911764]. [PubMed Central: PMC5599304].
- Shi T, McAllister DA, O'Brien KL, Simoes EAF, Madhi SA, Gessner BD, et al. Global, regional, and national disease burden estimates of acute lower respiratory infections due to respiratory syncytial virus in young children in 2015: A systematic review and modelling study. *Lancet.* 2017;**390**(10098):946–58. doi: 10.1016/S0140-6736(17)30938-8. [PubMed: 28689664]. [PubMed Central: PMC5592248].
- Caballero MT, Polack FP. Respiratory syncytial virus is an "opportunistic" killer. *Pediatr Pulmonol*. 2018;53(5):664-7. doi: 10.1002/ppul.23963. [PubMed: 29461021]. [PubMed Central: PMC5947624].
- Karron RA, Zar HJ. Determining the outcomes of interventions to prevent respiratory syncytial virus disease in children: What to measure? Lancet Respir Med. 2018;6(1):65–74. doi: 10.1016/S2213-2600(17)30303-X. [PubMed: 28865676].
- Zanone SM, Krause LK, Madhi SA, Bassat Q, Jha P, Simoes EA, et al. Challenges in estimating RSV-associated mortality rates. *Lancet Respir Med.* 2016;4(5):345–7. doi: 10.1016/S2213-2600(16)30042-X. [PubMed: 27304556].
- Carroll KN, Gebretsadik T, Escobar GJ, Wu P, Li SX, Walsh EM, et al. Respiratory syncytial virus immunoprophylaxis in high-risk infants and development of childhood asthma. J Allergy Clin Immunol. 2017;**139**(1):66–71 e3. doi: 10.1016/j.jaci.2016.01.055. [PubMed: 27212083]. [PubMed Central: PMC5074917].
- de Velde F, Mouton JW, de Winter BCM, van Gelder T, Koch BCP. Clinical applications of population pharmacokinetic models of antibiotics: Challenges and perspectives. *Pharmacol Res.* 2018;**134**:280–8. doi: 10.1016/j.phrs.2018.07.005. [PubMed: 30033398].
- 8. Shah JN, Chemaly RF. Management of RSV infections in adult recipients of hematopoietic stem cell transplantation. *Blood.*

2011;**117**(10):2755-63. doi: 10.1182/blood-2010-08-263400. [PubMed: 21139081].

- Pilie P, Werbel WA, Riddell J, Shu X, Schaubel D, Gregg KS. Adult patients with respiratory syncytial virus infection: Impact of solid organ and hematopoietic stem cell transplantation on outcomes. *Transpl Infect Dis.* 2015;17(4):551-7. doi: 10.1111/tid.12409. [PubMed: 26059180].
- Robinson JL, Grenier D, MacLusky I, Allen UD. Respiratory syncytial virus infections in pediatric transplant recipients: A Canadian Paediatric Surveillance Program study. *Pediatr Transplant*. 2015;19(6):659– 62. doi: 10.1111/petr.12553. [PubMed: 26152857].
- Griffiths C, Drews SJ, Marchant DJ. Respiratory syncytial virus: Infection, detection, and new options for prevention and treatment. *Clin Microbiol Rev.* 2017;**30**(1):277-319. doi: 10.1128/CMR.00010-16. [PubMed: 27903593]. [PubMed Central: PMC5217795].
- Wang M. Research progress of respiratory tract infection induced by respiratory syncytial virus in children. Int J Pediatr. 2018;45(7):501–4.
- Altman DG. Practical statistics for medical research. London: Chapman & Hall/CRC; 1991.
- Yun SG, Kim MY, Choi JM, Lee CK, Lim CS, Cho Y, et al. Comparison of three multiplex PCR assays for detection of respiratory viruses: Anyplex II RV16, AdvanSure RV, and Real-Q RV. J Clin Lab Anal. 2018;32(2). doi: 10.1002/jcla.22230. [PubMed: 28397965]. [PubMed Central: PMC5836940].
- Kanwar N, Hassan F, Nguyen A, Selvarangan R. Head-to-head comparison of the diagnostic accuracies of BD veritor system RSV and quidel(R) sofia(R) RSV FIA systems for respiratory syncytial virus (RSV) diagnosis. *J Clin Virol*. 2015;65:83–6. doi: 10.1016/j.jcv.2015.02.008. [PubMed: 25766995].
- Bouzas ML, Oliveira JR, Queiroz A, Fukutani KF, Barral A, Rector A, et al. Diagnostic accuracy of digital RNA quantification versus real-time PCR for the detection of respiratory syncytial virus in nasopharyngeal aspirates from children with acute respiratory infection. *J Clin Virol.* 2018;**106**:34–40. doi: 10.1016/j.jcv.2018.07.003. [PubMed: 30031351].