



Evaluation of Two Rapid Diagnostic *Clostridioides difficile* Infection Tests in a Chinese Hospital: A Real-world Analysis

Ge Huang¹, Yizheng Zhou¹, Tao Lv², Lisi Zheng², Yue Pei³, Yunbo Chen^{2,4,*} and Chengbin Li^{1,**}

¹Jingzhou Hospital Affiliated to Yangtze University, Jingzhou, China

²State Key Laboratory for Diagnosis and Treatment of Infectious Diseases, Collaborative Innovation Center for Diagnosis and Treatment of Infectious Diseases, The First Affiliated Hospital, School of Medicine, Zhejiang University, Hangzhou, China

³Clinical Laboratory, Women's Hospital, School of Medicine, Zhejiang University, Hangzhou, China

⁴Jinan Microecological Biomedicine Shandong Laboratory, Jinan, China

*Corresponding author: State Key Laboratory for Diagnosis and Treatment of Infectious Diseases, Collaborative Innovation Center for Diagnosis and Treatment of Infectious Diseases, The First Affiliated Hospital, School of Medicine, Zhejiang University, Hangzhou, China. Email: chenyunbo@aliyun.com

**Corresponding author: Jingzhou Hospital Affiliated to Yangtze University, Jingzhou, China. Email: yzchengbinli@sina.com

Received 2022 June 20; Revised 2022 October 22; Accepted 2022 October 25.

Abstract

Background: Accurate diagnosis is essential for optimal prevention and treatment of *Clostridioides difficile* infection (CDI), and various diagnostic methods must be evaluated.

Objectives: We aimed to evaluate and compare the performance of VIDAS *C. difficile*, C. DIFF QUIK CHEK COMPLETE (QCC), and toxigenic culture (TC) tests for diagnosing CDI and further determine the relationships between clinical factors and the toxin status of patients.

Methods: Stool samples were randomly selected for VIDAS or QCC testing according to the manufacturer's instructions between May 2017 and May 2021, and their performance was compared with that of TC. Clinical information was obtained from the hospital's electronic medical records.

Results: Among 10,897 samples tested, 6,435 and 4,462 samples were assigned for VIDAS and QCC tests, respectively. A total of 9.1% (996/10,897) of the samples were positive for TC. The sensitivity, specificity, positive predictive value, and negative predictive value were 36.6%, 98.6%, 72.1%, and 87.6% for VIDAS toxins A and B testing and 31.6%, 98.2%, 64.0%, and 87.8% for QCC toxin testing, respectively. Our results showed that the clinical data of the patients with positive and detectable toxins were not significantly different.

Conclusions: The VIDAS and QCC tests provide rapid screening assays for the laboratory diagnosis of CDI. However, a more specific test to detect free toxins is required to confirm the diagnosis for glutamate dehydrogenase (GDH)-positive and toxin-negative samples. The clinical characteristics and outcomes of this cohort were similar, regardless of the results of toxins A and B testing.

Keywords: *Clostridioides difficile*, Clinical Laboratory Diagnosis, Enzyme Immunoassay, Glutamate Dehydrogenase

1. Background

Clostridioides difficile is an anaerobically grown Gram-positive toxigenic clostridium that is the primary pathogen of antibiotic-associated diarrhea (1). The clinical manifestations of *C. difficile* infection (CDI) range from asymptomatic carriage to severe life-threatening toxic megacolon, sepsis, and death (2). The 2019 Antibiotic Resistance Threats report by the United States Centers for Disease Control and Prevention listed *C. difficile* as an urgent threat, as it caused more than 12,000 deaths per year in the United States (3). The incidence of *C. difficile* in developing countries was 8.5 cases per 10,000 patient days, which was approximately twice that in European countries. The situation is even worse considering the underdiagnosis of *C. difficile* in developing countries (4).

Due to the use of broad-spectrum antibiotics in recent years, the detection rate of toxin-producing *C. difficile* in China has increased from 8% to 14.8% (5-7).

The reference methods for diagnosing CDI are toxigenic culture (TC) and cell culture toxin neutralization assay (8). However, given the labor-intensive and time-consuming nature of the test, its use in clinical laboratories has been limited. Accurate diagnosis is essential for the optimal treatment and prevention of CDI, and a two-step test algorithm, which comprises a screening test with high sensitivity followed by a more specific test, has been recommended by some guidelines (9, 10). Direct fecal detection of *C. difficile* antigen glutamate dehydrogenase (GDH) and toxins antigens replaces traditional culture methods. Currently, the VIDAS GDH and toxin A/B

and C. DIFF QUIK CHEK COMPLETE (C. DIFF QCC; hereafter, QCC) tests are the common commercially available assays for detecting CDI in China. VIDAS GDH and toxins A and B assays are based on the enzyme-linked fluorescence immunoassay-related toxin detector enzyme-linked immunosorbent assay (EIA) method. QCC is a rapid diagnostic assay that combines the GDH antigen and toxins A and B detections.

2. Objectives

The awareness of CDI is increasing among Chinese clinicians. However, a methodological evaluation has not yet been established in China to help clinicians diagnose CDI accurately (11, 12). Therefore, we aimed to evaluate the performance of VIDAS and QCC tests relative to TC for the diagnosis of CDI and further assess the relationships between clinical factors and different toxin statuses of patients in this real-world study.

3. Methods

3.1. Sampling

This retrospective study was conducted at a 2200-bed tertiary care teaching hospital in China. Consecutive specimens from patients aged ≥ 18 years with clinically significant diarrhea (taking the shape of the receptacle or corresponding to Bristol Stool Type 5–7 plus ≥ 3 episodes of diarrhea within 24 h), excluding laxative-induced diarrheal stools, stools from patients who were repeatedly tested within seven days, and stools with a transmission time longer than 72 h, were sent to the clinical laboratory. Only liquid or unformed stools were processed using either the VIDAS GDH and Toxins A&B kit (bioMérieux, France) or C. DIFF QUIK CHEK COMPLETE kit (Techlab, Blacksburg, VA, USA) from May 2017 to May 2021. All samples were also tested using TC as the gold standard.

3.2. Definition

Patients suffering from diarrhea without other explanations and with the *C. difficile* toxin detected in their fecal samples were diagnosed with CDI, as defined previously (2). The severity score of CDI was based on the clinical practice guidelines of the Society for Healthcare Epidemiology of America and the Infectious Diseases Society of America (13). Severe CDI was characterized by a leukocyte count of $> 15,000$ cells/mL and a serum creatinine level of > 1.5 mg/dL.

3.3. Diagnostic Tests

Our laboratory has developed exclusion criteria for stool specimens. Therefore, stool samples that were formed or those that did not take the shape of the vessel upon inversion were excluded. The stool samples were tested for GDH and toxins using the VIDAS *C. difficile* GDH and toxins A&B kits (bioMérieux, France), respectively, as previously described (14). VIDAS toxins A&B testing was performed only if the VIDAS GDH test was positive. The QCC tests (Techlab Inc., Blacksburg, VA, USA) for both GDH (QCC-GDH) and toxins A/B (QCC-Toxins A&B) in a single cartridge were used to assess the stool samples according to the manufacturer's instructions (15). We added an equivalent volume of stool samples to a diluent and conjugated them in a tube (TechLab). Then, we transferred the mixture to the well of the device sample, followed by incubation for 15 min at room temperature. Next, we added the wash buffer and substrate (TechLab) to the reaction window and read the result after 10 min. The GDH antigen and/or toxins were reported as positive if a corresponding band was seen on the device window.

3.4. Gold Standard

Unformed stools were processed immediately or, if logistically impractical, stored at 4°C for 24 h until processing. All samples were inoculated anaerobically on cycloserine-cefoxitin fructose agar (Oxford, UK) to isolate *C. difficile* using alcohol shock after VIDAS or QCC testing. Colonies of *C. difficile* were confirmed by both matrix-assisted laser desorption ionization-time of flight mass spectrometry (MALDI/TOF MS, Bruker Corporation, USA). The *tcdA*, *tcdB*, and *cdtA/B* genes in the isolated *C. difficile* strains were detected by PCR, as described elsewhere (6). If the *C. difficile* isolates were positive on PCR for any of the *tcdA*, *tcdB*, and *cdtA/B* genes, the specimen was graded as TC-positive. The clinicians were blinded to this result.

3.5. Clinical Data

We performed a retrospective cohort study evaluating the TC-positive stool samples and clinical characteristics of patients with positive VIDAS and QCC test results for both GDH and toxins A and B (toxin-positive group) and those with GDH-positive and toxin-negative results (toxin-negative/TC-positive group). The clinical information obtained included age, sex, oral antibiotic (vancomycin or metronidazole) use before sample submission, basic diseases, routine clinical blood tests on the day of diarrhea onset or three days before and after, and clinical prognostic results (including ICU admissions and death). All data were obtained from the electronic medical records of the hospital.

3.6. Statistical Analysis

The sensitivity and specificity of the assays for GDH and toxins A&B were calculated against the results of TC as the gold standard. Cohen's kappa was computed to evaluate the inter-assay agreement between the VIDAS and QCC tests (agreement: < 0.4, poor; 0.4 - 0.75, fair to good; > 0.75, excellent). The chi-squared or Fisher's exact test was used to compare the categorical variables, and the Mann-Whitney U test and Wilcoxon rank sum test were used to compare the continuous variables. The specificity, sensitivity, positive predictive value (PPV), and negative predictive value (NPV) were also computed. The confidence intervals (CIs) were calculated at the 95% level, and the kappa value for the test performance evaluation in this comparative study was determined using SPSS software (SPSS 25.0, IBM Corporation, Armonk, NY, USA). For all analyses, two-tailed P-values of < 0.05 denoted statistical significance.

4. Results

4.1. General Performance

We evaluated 12,268 stool specimens submitted for *C. difficile* detection between May 2017 and May 2021. After the exclusion of duplicate specimens and unformed stools, 10,897 stool specimens met the inclusion criteria for this study. All samples were randomly selected for VIDAS (6,435) or QCC (4,462) testing. Of 10,897 samples, 996 (9.1%) tested positive for TC. For VIDAS detection, 297 (4.6%) patients tested positive for GDH and toxins A and B; 933 (14.5%) patients were GDH-positive and toxins A and B-negative, and a total of 5,205 (80.9%) patients tested negative for GDH and toxins A and B. On the other hand, 289 (31.0%) of 933 patients and 81 (1.6%) of 5,205 patients tested positive for TC (Figure 1). Of the 4,462 samples, 203 (4.5%) were GDH- and toxin-positive, and 3,613 (81.0%) were GDH-positive and toxin-negative. Of the 646 GDH-positive/toxin-negative samples tested using TC, 224 (34.7%) were positive. Of the 3,613 GDH-positive/toxin-negative samples, 58 (1.8%) were positive, while the remaining were negative on TC (Figure 2).

4.2. Sensitivity, Specificity, Positive Predictive Value, and Negative Predictive Value of Glutamate Dehydrogenase Detection Assay

The sensitivity and specificity of the VIDAS GDH test were 86.1% and 98.4%, while those for the QCC GDH test were 85.9% and 98.4%, respectively (Table 1). Both tests showed lower PPVs than TC, which were 40.9% and 41.7% for VIDAS GDH and QCC GDH, respectively.

4.3. Sensitivity, Specificity, Positive Predictive Value, and Negative Predictive Value of Toxin Detection Assays

Using TC as a reference method, the sensitivities and specificities of the tests were as follows: 36.6% sensitivity and 98.6% specificity for VIDAS toxins A&B, and 31.6% sensitivity and 98.2% specificity for the QCC toxin test (Table 1). Both tests showed moderate PPVs, which were 72.1% and 64.0% for VIDAS toxins, A&B, and QCC toxins, respectively. The NPVs of the VIDAS toxins A&B and QCC toxin tests were similar (95.3% and 95.1%, respectively).

4.4. Discrepant Results of Toxigenic Culture

Eighty-three patient specimens tested positive on the VIDAS GDH and toxins A&B assays and negative on TC, while 73 specimens tested positive on the QCC GDH and toxin assays but negative on TC.

4.5. Clinical Data

A total of 399 patients with diarrhea who tested positive for *C. difficile* on TC and had complete clinical information were included in the study, which included 241 patients tested using VIDAS and 158 patients tested using QCC. In the VIDAS group, 90 patients were classified as the GDH-positive/toxins A and B-positive group and 151 as the GDH-positive/toxins A and B-negative group. There were no statistically significant differences in other demographic variables, including age and sex (Appendix 1 in Supplementary File). There was no significant difference in disease distribution between the two groups. The two groups did not have significantly different clinical outcomes (mortality) and unfavorable prognoses (ICU admissions and severe CDI) ($P > 0.05$). Twenty-four patients, including nine in the GDH-positive/toxins A&B-negative group and 15 in the GDH-positive/toxins A and B-negative group, were treated with oral vancomycin or metronidazole before testing.

In the QCC group, 56 patients belonged to the GDH-positive/toxins-positive group, and 102 belonged to the GDH-positive/toxins-negative group. The data on age, sex, predisposing factors, ICU admission, creatinine values, and mortality were similar in the two groups (Appendix 2 in Supplementary File). Similarly, we observed two GDH-positive/toxin-negative and 11 GDH-positive/toxin-negative stool samples submitted after receiving clinically empirical treatment at the time of diarrhea onset; however, these differences were not statistically significant ($P = 0.140$).

5. Discussion

With the incidence and severity increasing annually (16), accurate and rapid diagnosis of CDI has become particularly important. In this real-world study, we found that




Figure 1. Flow chart of VIDAS detection of toxigenic *Clostridioides difficile*. GDH, VIDAS-GDH; TC, toxigenic culture.




Figure 2. Flow chart of QCC detection of toxigenic *Clostridioides difficile*. GDH*, C. DIFF QCC antigen GDH; toxins, C. DIFF QCC-toxin; TC, toxigenic culture.

Table 1. Sensitivity and Specificity of Algorithms Compared with the Reference Method ^a

Diagnostic Method	No. of Specimens				%Sensitivity (95%CI)	%Specificity (95%CI)	%PPV (95%CI)	%NPV (95%CI)	Kappa
	TP	FP	TN	FN					
GDH+A&B (VIDAS)	214	83	5768	370	36.6 (32.7 - 40.6)	98.6 (98.3 - 98.9)	72.1 (66.9 - 77.2)	94.0 (93.4 - 94.6)	0.452
GDH ⁺ +toxins (QCC)	130	73	3977	282	31.6 (27.0 - 36.1)	98.2 (97.8 - 98.6)	64.0 (57.4 - 70.7)	93.4 (92.6 - 94.1)	0.39

^a GDH, VIDAS-GDH; A&B, VIDAS A&B; GDH⁺, QCC antigen GDH; toxins, C. DIFF QCC-toxin; TP, true positive; FP, false positive; TN, true negative; PPV, positive predictive value; NPV, negative predictive value; CI, confidence interval.

GDH detection using the VIDAS and QCC kits was highly sensitive and consistent with the manufacturer's instructions. On the other hand, the sensitivity of toxin detection was significantly lower than that described in the kit instructions or other reports (12, 17). Furthermore, we found that the clinical characteristics and outcomes of this cohort were similar regardless of toxin status, emphasizing the clinical evaluation of patient status (18). Therefore, a negative GDH can exclude the possibility of CDI, while in patients with inconsistent GDH and toxin results, further tests, such as culture, are needed to identify the toxin.

In the present study, the QCC and VIDAS tests showed favorable performance for detecting GDH with similar sensitivity and specificity to the previously reported data (87.5 - 95.8% and 82.9 - 97.5%, respectively) (19, 20). In areas with low CDI prevalence, a high NPV can reliably rule out the diagnosis of CDI (21). However, the sensitivities of the VIDAS and QCC tests for the detection of toxins were 36% and 31.6%, respectively, which were lower than the range reported in previous studies (sensitivities, 51 - 60%; specificities, 93.4 - 97.4%) (17, 22). This may be attributed to the different reference methods employed (23). Another reason may be the different prevalent strains in different regions or settings. For example, toxin-hyperproducing NAP1 strains are more often associated with toxin-positive results (24). Given the low sensitivity of EIA, additional tests such as the Nucleic Acid Amplification Test (NAAT) should be performed when CDI is suspected, but the EIA is negative (25). With this, a misdiagnosis of CDI may be avoided (26).

There is still controversy about the need to treat patients with evidence of toxigenic *C. difficile* but negative toxins A&B based on EIA, as they may have CDI with undetectable toxin levels, false-negative toxins A&B EIA, or asymptomatic colonization (27). We compared the outcomes of patients with toxin EIA-negative, TC-positive, and toxin EIA-positive/TC-positive results to assess their clinical characteristics. The results showed no difference in clinical characteristics and outcomes, consistent with previous reports (28, 29). However, research in the US reported that patients with positive NAAT-toxins had more complications and a higher blood leucocyte count than patients who were only positive for NAAT (30). This may be due to the prevalent strains in our setting; non-binary-toxin-producing *C. difficile*, which causes less severe CDI, was the

major toxin type in our setting, while NAP1/BI/027 was the common type in the US (31, 32). In our study, severe CDI occurred in 15.9% of the VIDAS toxin-negative and 7.6% of QCC toxin-negative patient groups, respectively. Toxin-EIA is unreliable in differentiating CDI (33), and we suggest that further testing, such as culture, is needed to detect toxins in stools where GDH is inconsistent with toxins.

Our study has some limitations. First, it was carried out in a single setting, and the findings may not apply to other centers with different demographics. Second, we used TC as the gold standard and detected toxin genes using PCR. Bacteria that harbor a toxin gene do not always produce the toxin (34). This may have reduced the sensitivity of toxin detection. Lastly, although our laboratory developed criteria to exclude formed stool specimens, however, there was no confirmation of clinical symptoms, which may result in the colonization of *C. difficile* rather than infection in some patients.

5.1. Conclusions

The VIDAS and QCC tests have similar specificities for diagnosing CDI, and both methods have good NPV for the rapid screening of negative patients. However, they have poor sensitivity for diagnosing CDI. More sensitive assays are needed for GDH-positive and toxin-negative samples to detect toxins for diagnostic confirmation.

Supplementary Material

Supplementary material(s) is available [here](#) [To read supplementary materials, please refer to the journal website and open PDF/HTML].

Footnotes

Authors' Contribution: Yunbo Chen and Chengbin Li designed the study. Ge Huang and Yizheng Zhou conducted the study. Tao Lv, Lisi Zheng, and Yue Pei completed the experiments. Ge Huang analyzed the data and wrote the paper. Yunbo Chen and Chengbin Li reviewed and edited the paper. All authors approved the final version of the article.

Conflict of Interests: The authors declare no conflicts of interest.

Ethical Approval: This study was approved by the hospital research ethics committee with an informed consent waiver (No.2020-651).

Funding/Support: This research was funded by the Jinan Microecology Laboratory (JNL-2022037C).

Informed Consent: The data are anonymous, and the requirement for informed consent was waived.

References

- Kocielek LK, Gerding DN. Breakthroughs in the treatment and prevention of Clostridium difficile infection. *Nat Rev Gastroenterol Hepatol*. 2016;13(3):150–60. [PubMed ID: 26860266]. <https://doi.org/10.1038/nrgastro.2015.220>.
- McDonald LC, Gerding DN, Johnson S, Bakken JS, Carroll KC, Coffin SE, et al. Clinical practice guidelines for Clostridium difficile infection in adults and children: 2017 update by the Infectious Diseases Society of America (IDSA) and Society for Healthcare Epidemiology of America (SHEA). *Clin Infect Dis*. 2018;66(7):e1–48. [PubMed ID: 29462280]. [PubMed Central ID: PMC6018983]. <https://doi.org/10.1093/cid/cix1085>.
- Centers for Disease Control and Prevention. Antibiotic resistance threats in the United States, 2019. 2019. Available from: <https://stacks.cdc.gov/view/cdc/82532>.
- Zyoud SH. Global research on Clostridium difficile-associated diarrhoea: A visualized study. *World J Gastroenterol*. 2022;28(28):3720–31. [PubMed ID: 36161039]. [PubMed Central ID: PMC9372798]. <https://doi.org/10.3748/wjg.v28.i28.3720>.
- Li C, Li Y, Huai Y, Liu S, Meng X, Duan J, et al. Incidence and outbreak of healthcare-onset healthcare-associated clostridioides difficile infections among intensive care patients in a large teaching hospital in China. *Front Microbiol*. 2018;9:566. [PubMed ID: 29636739]. [PubMed Central ID: PMC5880889]. <https://doi.org/10.3389/fmicb.2018.00566>.
- Chen YB, Gu SL, Wei ZQ, Shen P, Kong HS, Yang Q, et al. Molecular epidemiology of Clostridium difficile in a tertiary hospital of China. *J Med Microbiol*. 2014;63(Pt 4):562–9. [PubMed ID: 24344206]. <https://doi.org/10.1099/jmm.0.068668-0>.
- Meng X, Huang X, Peng Z, Wang Y, Liu S, Zeng C, et al. Antibiotic resistances and molecular characteristics of Clostridioides difficile in ICUs in a teaching hospital from central south China. *Front Med (Lausanne)*. 2021;8:745383. [PubMed ID: 34938744]. [PubMed Central ID: PMC8685378]. <https://doi.org/10.3389/fmed.2021.745383>.
- Crobach MJT, Baktash A, Duszynski N, Kuijper EJ. Diagnostic guidance for C. difficile infections. *Adv Exp Med Biol*. 2018;1050:27–44. [PubMed ID: 29383662]. https://doi.org/10.1007/978-3-319-72799-8_3.
- Johnson S, Lavergne V, Skinner AM, Gonzales-Luna AJ, Garey KW, Kelly CP, et al. Clinical practice guideline by the Infectious Diseases Society of America (IDSA) and Society for Healthcare Epidemiology of America (SHEA): 2021 focused update guidelines on management of clostridioides difficile infection in adults. *Clin Infect Dis*. 2021;73(5):e1029–44. [PubMed ID: 34164674]. <https://doi.org/10.1093/cid/ciab549>.
- Kelly CR, Fischer M, Allegretti JR, LaPlante K, Stewart DB, Limketkai BN, et al. ACG clinical guidelines: Prevention, diagnosis, and treatment of Clostridioides difficile infections. *Am J Gastroenterol*. 2021;116(6):1124–47. [PubMed ID: 34003176]. <https://doi.org/10.14309/ajg.0000000000001278>.
- Moon HW, Kim HN, Hur M, Shim HS, Kim H, Yun YM. Comparison of diagnostic algorithms for detecting toxigenic Clostridium difficile in routine practice at a tertiary referral hospital in Korea. *PLoS One*. 2016;11(8):e0161139. [PubMed ID: 27532104]. [PubMed Central ID: PMC4988646]. <https://doi.org/10.1371/journal.pone.0161139>.
- Savolainen R, Koskinen JM, Mentula S, Koskinen JO, Kaukoranta SS. Prospective evaluation of the mapiPOC test for detection of Clostridioides difficile glutamate dehydrogenase and toxins A/B. *J Clin Microbiol*. 2020;58(4):e01872–19. [PubMed ID: 31941691]. [PubMed Central ID: PMC7098770]. <https://doi.org/10.1128/JCM.01872-19>.
- Cohen SH, Gerding DN, Johnson S, Kelly CP, Loo VG, McDonald LC, et al. Clinical practice guidelines for Clostridium difficile infection in adults: 2010 update by the society for healthcare epidemiology of America (SHEA) and the infectious diseases society of America (IDSA). *Infect Control Hosp Epidemiol*. 2010;31(5):431–55. [PubMed ID: 20307919]. <https://doi.org/10.1086/651706>.
- Shin BM, Lee EJ, Moon JW, Lee SY. Evaluation of the VIDAS glutamate dehydrogenase assay for the detection of Clostridium difficile. *Anaerobe*. 2016;40:68–72. [PubMed ID: 27282799]. <https://doi.org/10.1016/j.anaerobe.2016.06.001>.
- Quinn CD, Sefters SE, Babiker W, He Y, Alcubasa R, Stratton CW, et al. C. Diff Quik Chek complete enzyme immunoassay provides a reliable first-line method for detection of Clostridium difficile in stool specimens. *J Clin Microbiol*. 2010;48(2):603–5. [PubMed ID: 19955275]. [PubMed Central ID: PMC2815622]. <https://doi.org/10.1128/JCM.01614-09>.
- Park H, Laffin MR, Jovel J, Millan B, Hyun JE, Hotte N, et al. The success of fecal microbial transplantation in Clostridium difficile infection correlates with bacteriophage relative abundance in the donor: a retrospective cohort study. *Gut Microbes*. 2019;10(6):676–87. [PubMed ID: 30866714]. [PubMed Central ID: PMC6867182]. <https://doi.org/10.1080/19490976.2019.1586037>.
- Yoo IY, Song DJ, Huh HJ, Lee NY. Simultaneous detection of Clostridioides difficile glutamate dehydrogenase and toxin A/B: Comparison of the C. DIFF QUIK CHEK COMPLETE and RIDASCREEN assays. *Ann Lab Med*. 2019;39(2):214–7. [PubMed ID: 30430786]. [PubMed Central ID: PMC6240517]. <https://doi.org/10.3343/alm.2019.39.2.214>.
- Dubberke ER, Han Z, Bobo L, Hink T, Lawrence B, Copper S, et al. Impact of clinical symptoms on interpretation of diagnostic assays for Clostridium difficile infections. *J Clin Microbiol*. 2011;49(8):2887–93. [PubMed ID: 21697328]. [PubMed Central ID: PMC3147743]. <https://doi.org/10.1128/JCM.00891-11>.
- Carson KC, Riley TV. Comparison of the Vidas C. difficile and Quik Chek-60 glutamate dehydrogenase assays for the detection of Clostridium difficile in faecal samples. *Pathology*. 2016;48(5):506–8. [PubMed ID: 27311869]. <https://doi.org/10.1016/j.pathol.2016.04.007>.
- Kim H, Kim WH, Kim M, Jeong SH, Lee K. Evaluation of a rapid membrane enzyme immunoassay for the simultaneous detection of glutamate dehydrogenase and toxin for the diagnosis of Clostridium difficile infection. *Ann Lab Med*. 2014;34(3):235–9. [PubMed ID: 24790912]. [PubMed Central ID: PMC3999323]. <https://doi.org/10.3343/alm.2014.34.3.235>.
- Chen YB, Gu SL, Shen P, Lv T, Fang YH, Tang LL, et al. Molecular epidemiology and antimicrobial susceptibility of Clostridium difficile isolated from hospitals during a 4-year period in China. *J Med Microbiol*. 2018;67(1):52–9. [PubMed ID: 29160203]. <https://doi.org/10.1099/jmm.0.000646>.
- Kim N, Lee SY, Park J, Lee J. Comparative evaluation of three immunoassays for the simultaneous detection of Clostridioides difficile glutamate dehydrogenase and toxin A/B. *Microorganisms*. 2022;10(5):947. [PubMed ID: 35630390]. [PubMed Central ID: PMC9145049]. <https://doi.org/10.3390/microorganisms10050947>.
- Bagdasarjan N, Rao K, Malani PN. Diagnosis and treatment of Clostridium difficile in adults: a systematic review. *JAMA*. 2015;313(4):398–408. [PubMed ID: 25626036]. [PubMed Central ID: PMC6561347]. <https://doi.org/10.1001/jama.2014.17103>.
- Gomez EJ, Montgomery S, Alby K, Robinson DP, Roundtree SS, Blecker-Shelly D, et al. Poor yield of Clostridium difficile testing algorithms using glutamate dehydrogenase antigen and C. difficile toxin enzyme immunoassays in a pediatric population with

- declining prevalence of clostridium difficile strain BI/NAP1/027. *Diagn Microbiol Infect Dis*. 2018;**91**(3):229–32. [PubMed ID: [29567127](#)]. <https://doi.org/10.1016/j.diagmicrobio.2018.02.016>.
25. Yoldas O, Altindis M, Cufali D, Asik G, Kesli R. A diagnostic algorithm for the detection of Clostridium difficile-associated diarrhea. *Balkan Med J*. 2016;**33**(1):80–6. [PubMed ID: [26966622](#)]. [PubMed Central ID: [PMC4767316](#)]. <https://doi.org/10.5152/balkanmedj.2015.15159>.
 26. Guery B, Galperine T, Barbut F. Clostridioides difficile: Diagnosis and treatments. *BMJ*. 2019;**366**:l4609. [PubMed ID: [31431428](#)]. <https://doi.org/10.1136/bmj.l4609>.
 27. Fang FC, Polage CR, Wilcox MH. Point-counterpoint: What is the optimal approach for detection of Clostridium difficile infection? *J Clin Microbiol*. 2017;**55**(3):670–80. [PubMed ID: [28077697](#)]. [PubMed Central ID: [PMC5328433](#)]. <https://doi.org/10.1128/JCM.02463-16>.
 28. Miller R, Morillas JA, Brizendine KD, Fraser TG. Predictors of Clostridioides difficile infection-related complications and treatment patterns among nucleic acid amplification test-positive/toxin enzyme immunoassay-negative patients. *J Clin Microbiol*. 2020;**58**(3):e01764–19. [PubMed ID: [31896665](#)]. [PubMed Central ID: [PMC7041579](#)]. <https://doi.org/10.1128/JCM.01764-19>.
 29. Bernard R, Hammami MB, Arnold FW, McGrath B, Patel A, Wuerth B, et al. Clostridioides difficile toxin is infrequently detected in inflammatory bowel disease and does not associate with clinical outcomes. *Gut Pathog*. 2022;**14**(1):36. [PubMed ID: [36042530](#)]. [PubMed Central ID: [PMC9426007](#)]. <https://doi.org/10.1186/s13099-022-00511-2>.
 30. Polage CR, Gyorke CE, Kennedy MA, Leslie JL, Chin DL, Wang S, et al. Overdiagnosis of Clostridium difficile infection in the molecular test era. *JAMA Intern Med*. 2015;**175**(11):1792–801. [PubMed ID: [26348734](#)]. [PubMed Central ID: [PMC4948649](#)]. <https://doi.org/10.1001/jamainternmed.2015.4114>.
 31. Xu Q, Chen Y, Gu S, Lv T, Zheng B, Shen P, et al. Hospital-acquired Clostridium difficile infection in Mainland China: A seven-year (2009-2016) retrospective study in a large university hospital. *Sci Rep*. 2017;**7**(1):9645. [PubMed ID: [28852010](#)]. [PubMed Central ID: [PMC5575102](#)]. <https://doi.org/10.1038/s41598-017-09961-0>.
 32. Tickler IA, Goering RV, Whitmore JD, Lynn AN, Persing DH, Tenover FC, et al. Strain types and antimicrobial resistance patterns of Clostridium difficile isolates from the United States, 2011 to 2013. *Antimicrob Agents Chemother*. 2014;**58**(7):4214–8. [PubMed ID: [24752264](#)]. [PubMed Central ID: [PMC4068552](#)]. <https://doi.org/10.1128/AAC.02775-13>.
 33. Parnell JM, Fazili I, Bloch SC, Lacy DB, Garcia-Lopez VA, Bernard R, et al. Two-step testing for clostridioides difficile is inadequate in differentiating infection from colonization in children. *J Pediatr Gastroenterol Nutr*. 2021;**72**(3):378–83. [PubMed ID: [32925555](#)]. [PubMed Central ID: [PMC7870537](#)]. <https://doi.org/10.1097/MPG.0000000000002944>.
 34. Morris KA, Davies K, Wilcox MH. Impact of Clostridium difficile toxin gene PCR result on decisions to de-isolate patients: Do the ends justify the means? *J Infect Prev*. 2018;**19**(3):138–40. [PubMed ID: [29796096](#)]. [PubMed Central ID: [PMC5956705](#)]. <https://doi.org/10.1177/1757177418755309>.